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## RT41

### GB STEAM AUTOCLAVE

Pre-set, non-modifiable cycles: all cycles are pre-set at the factory and guarantee perfect and safe sterilization of treated materials. There is no way the operator can manipulate or alter the cycles, and this safeguards against risks deriving from human error-Pre-set temperatures: the temperatures set for RT41 autoclave cycles are 134°C and 121°C. The 134°C cycles are normally used for the sterilization of solid and ferrous materials. 121°C cycles are normally used for sterilization of delicate materials-Operating cycles: RT41 is equipped to perform 2 sterilization cycles: 1 sterilization cycle at 134°C, 1 sterilization cycle at 121°C-Fully automatic preheating, sterilization and drying phases-Closure system with operator safety pin. Thermodynamic system for expelling "cold" air pockets-Chamber made entirely of stainless steel, with band-type heating element to guarantee even heating-4 aluminium trays.

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Loading classification / Classification des charges Type N

Weight / poids kg 35  
 Power consumption / Puissance W 1500  
 Voltage / voltage 230 V - 50/60 hz

Class Classe	N	16 Litres	Litres
1500 W	W	2 Cycles	Cycles

## RT92

### GB FRACTIONAL VACUUM STEAM AUTOCLAVE

Treated circular chamber made of AISI 304 stainless steel-Maximum load per cycle 5 kg-3 safety devices on the door: electromagnetic, electrical, dynamic-Controls card installed in the board-Graphic display-Automatic barometric alignment system-Pump-operated water filling system with automatic stop when the tank is full-Door seal to ensure maximum duration-Rear board with prefitings and automatic water filling and discharge-Forced ventilation system incorporated in the equipment-Protected tray holder to safeguard against contact between load and chamber-Fractional vacuum-RT92 is equipped with a fractional vacuum system. This system alternates negative pressure (vacuum) phases with positive pressure phases. The system is activated at the beginning of all sterilisation cycles. This method is the most technologically advanced system to produce vacuum and ensures the complete elimination of air pockets in all types of materials-Sterilization cycles-Pre-set temperatures: the preset temperatures used in our autoclaves' cycles are 134°C and 121°C. 134°C cycles are normally used for sterilization of solid and ferrous materials. 121°C cycles are normally used for sterilization of thermoplastic and delicate materials. In sterilizing materials, always follow the recommendations of the device manufacturer-Ease of use: choosing and selecting a sterilization cycle on RT92 autoclaves is extremely simple and intuitive. One of the buttons makes it possible to scroll and view data about all sterilization cycles on the display, thus allowing you to find the cycle that suits your materials best, then you will only have to push another button to automatically start the whole sterilization cycle without having to worry about anything. In case of any error or warning a series of messages and LEDs will attract your attention to the problem-Fully automatic cycles: all cycles are fully automatic, with all steps managed by the leading-edge computer board with a powerful microprocessor. Each step (vacuum, preheating, exposure time, drying) is automatically monitored, assisted and handled. There is no need for human intervention in cycle management-Wide choice of operation cycles: RT92 is equipped with: 1 sterilization cycle temperature 134°C for free or packed materials, 1 sterilization cycle temperature 121°C for free or packed materials, 1 "flash" cycle for quick sterilization cycles, 1 "prion" cycle especially designed for Creutzfeldt & Jacobs disease-Cycles without operator intervention and with auto-shutdown. All sterilization cycles have an automatic function that, if the operator is not by the machine when the cycle ends, will trigger an anti-condensation process and auto-shutdown operation. This system makes it possible to economize and increase safety of the device-Test cycles: our autoclaves are equipped with test cycles for periodic control of the device by means of physical and functional tests: Helix test to evaluate vacuum capability in porous bodies, Bowie & Dick test to evaluate vacuum capability in hollow bodies, Vacuum test to evaluate vacuum retaining capability.

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## DRY HEAT STERILIZING UNIT SMALL

### GB DRY HEAT STERILIZING UNIT SMALL

Outer case made of stainless steel. Inner structure, grid, removable perforated tray, hinged lid completely made of stainless steel. Inner structure and lid with nontoxic silicon rubber gasket, insulated by high grade glass wool. Stainless steel armoured heating element. Timer from 0 to 120 minutes, which disconnects the device at the end of the programmed time. Probe thermometer 0° - 200°C. Probe thermo-regulator 0° - 200°C. Safety thermostat. Red pilot light shows the line insertion. Green pilot light points out the resistance operation. 230V - 50/60 Hz.

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## DRY HEAT ELECTRIC STERILIZING UNIT

### GB DRY HEAT ELECTRIC STERILIZING UNIT

Outer case made of stainless steel. Inner structure, two adjustable inner shelves, hinged door completely made of stainless steel. Inner structure and door with nontoxic silicon rubber gasket, insulated by high grade glass wool. Stainless steel armoured heating element. Lockable with key. Adjustable aeration tunnels, placed on the bottom of the right side and on the top of the left side for hot air adjustment and circulation, which is useful for a quick cooling of the device. The electric plant is insulated and protected by silicon sheaths. Timer from 0 to 120 minutes, which disconnects the device at the end of the programmed time. Probe thermometer 0° - 200°C. Probe thermo-regulator 0° - 200°C. Safety thermostat. Double control pilot light (red and green) 230V - 50/60 Hz.

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CAPACITY / CAPACITÉ  
 A3213400 lt. 20  
 A3214535 lt. 60  
 A3215670 lt. 120

## DRY HEAT ELECTRIC STERILIZING UNIT FORCED HOT AIR

### GB DRY HEAT ELECTRIC STERILIZING UNIT FORCED HOT AIR

Outer case and front completely made of stainless steel, included the drummed door strongly isolated with non-conducting material. Inner chamber entirely made of stainless steel, strongly strengthened, resistant to the heat and the corrosion at hot temperature. Inner stainless steel shelves adjustable in height. Door with lock and key. To guarantee the uniformity of the temperature in the heating room and to allow the circulation of the air, the device is provided with a ventilation motor that results always in operation, both during the cycle of heating and during the pauses of cooling. The wiring is isolated and protected by silicon gaskets. Timer 0-120 minutes that switch off the oven at the end of the planned time. Probe thermometer from 0° to 200°C. Probe thermoregulator from 0° to 200°C. Safety thermostat. Ventilator for forced ventilation. Double check up light bulb. 230V-50/60 Hz.

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## DRY HEAT ELECTRONIC STERILIZING UNIT

### GB DRY HEAT ELECTRONIC STERILIZING UNIT

Outer case made of stainless steel. Inner structure, two adjustable inner shelves and hinged door completely made of stainless steel. Inner structure and door with nontoxic silicon rubber gasket, insulated by high grade glass wool. Stainless steel armoured heating element. Lockable by key. Adjustable aeration tunnels, placed on the bottom of the right side and on the top of the left side for hot air adjustment and circulation, which is useful for a quick cooling of the device. Electronic running C.P.S. Use modality: continuous (stand by) or timed. Timer from 0 to 250 minutes. Adjustable temperature from 0 to 200°C. Safety thermostat. Optimum heat distribution. Automatic indication for possible problems of feeding (in case of mains tension lack, the program will start again from the beginning). 230V – 50/60 Hz.

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### GB STERILIZATION REELS AND SELF-SEALING BAGS

These devices are designed for hospital, medical, dental, cosmetic, veterinary uses and, in general, for all the personnel that carries out steam or gas sterilization processes. Sterilization reels and bags are made of composite plastic material (polyester-polypropylene) and special medical paper, moisture resistant and with high anti-bacteria protection. The closure system of the bags prevent the forming of air bubbles and ensures a perfect seal, preventing fine dust from entering the bag and therefore causing unsatisfactory sterility. The inks used are water based. Rolls feature a toning indicator every 10 cm while each bag has one toning indicator. The indicators are to be used for steam or gas toning.

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#### SELF SEALING BAGS / ENVELOPPES AUTOSCELLANTES

SSP09025/1 90 mm x 250 mm (box of 200 pcs)  
 SSP14026/1 140 mm x 260 mm (box of 200 pcs)  
 SSP19040/1 190 mm x 330 mm (box of 200 pcs)  
 SSP30050/1 300 mm x 450 mm (box of 200 pcs)

## RT20

### GB THERMOSEALING MACHINE

RT20 is a new generation sealer that has been created with the aim of reducing envelope packaging time, monitoring the quality of the welding, and is able to seal all types of paper for sterilization, including folded paper-RT20 assures: sealing band of 12 mm, sealing bar 32 cm, great easiness of paper introduction, active electronic control over the couple of resistors, visibility of the area of cutting, signal of possible bad sealing, introducible in any sterilizing piece of furniture.

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## RT30

### GB ROLL SEALING MACHINE, WITH AUTOMATIC FEED SYSTEM WITH PHOTOCCELL

RT30 is a roll sealing machine, with automatic feed system with photocell. This is the maximum expression of technology applied to medical sealing, and is used in medical centers, dentists' surgeries, hospitals and clinics, ensuring very high output, maximum operating speed and no routine maintenance, and uses pre-cut bags sealed on 3 sides-RT30 is complete with integrated printer that, during sealing, automatically stamps the packaging date and expiry date, as well as the required regulatory symbols, on the bag-This very important solution allows you to always have control over the real sterility maintenance state of the instruments-On the side of the equipment there is also an electronic control and management card that allows the parameters of use to be adjusted-The use of this machine is simple and rational-All the operator has to do is insert the bag with the instrument on the sealing machine support and the photocell automatically and autonomously feeds the bag, seals it and stamps the packaging data-This action can be continuative and repetitive to ensure very high output and considerable reduction of time.

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